

# VENOSAN® AES



## *Thrombosis and embolism prevention* *Anti-embolism stockings*

The anti-embolism stockings are very easy to put on and offer a high level of comfort in wear due to their breathability.

### PERFECT FIT

Anatomical shaping and comprehensive choice of sizes.

### COLOUR-CODED

Allows easy and reliable size recognition after washing.



### SAFE

Because of the top inspection window.

# VENOSAN® AES *Thrombosis prevention and Anti-embolism stockings*



<b>SI</b>	<b>O</b>	Pressure at the ankle	<b>Material:</b>	<b>CE</b>
C A: S		<b>18 mmHg</b>	90% Nylon   Polyamide, 10% Spandex (Latex Free)	

Article & Swatch	Ankle Circumference	Calf Circumference	Thigh Circumference	Length	Size	Size colour code (toe)	Length colour code (top)	Art. No
<b>AD</b> Below knee 	18–20 cm	to 32 cm		to 40.5 cm	S regular	yellow	natural	SG57000
				over 40.5 cm	S long	yellow	blue	SG57001
	20–22 cm	32–37 cm		to 43 cm	M regular	natural	natural	SG57002
				over 43 cm	M long	natural	blue	SG57003
	23–25 cm	37–42 cm		to 45.5 cm	L regular	blue	natural	SG57004
				over 45.5 cm	L long	blue	blue	SG57005
<b>SG</b> White 	25–28 cm	42–47 cm		to 45.5 cm	XL regular	green	natural	SG57006
				over 45.5 cm	XL long	green	blue	SG57007
	28–30 cm	47–52 cm		to 45.5 cm	XXL regular	purple	natural	SG57008
				over 45.5 cm	XXL long	purple	blue	SG57009

**leef.**  
reverse how you feel

**Customer Service | Solutions & Sales**

1300 005 333 | [www.leef.com.au](http://www.leef.com.au) | [info@leef.com.au](mailto:info@leef.com.au)  
652 Glen Huntly Road, Caulfield South, Victoria 3162